

# 2510 Purging unit



Spülbecher 2500-0002-60

## Description

Automatically rotating purging unit.


A cup of delrin is inserted into a magnetic holder in the device and can be easily exchanged.

If the dosing needle is stripped on the edge of the cup, the container rotates by 3.6 °.  
A signal "container full" or "container not available" is output.

Space-saving unit of purging system and needle sensor.

# 2510 Purging unit

Dimensional drawing (Dimensions in mm)



Drive	Stepper motor DC 24V
Strip-off step	3,6°
Wearing part purging cup	Art. 2500-0002-60 / Delrin, 100 ml / Ø 70 mm
Dimensions	160 x 135 x 126 mm
Weight	1 kg